

imp De Mayo V²™ Knee Positioner



- ◆ Lighter Base Plate.
 - ◆ *New Carriage* with improved holding power.
- ◆ *New side Rail Clamp and Base Plate Bar* improving base plate security.



Illustrated with the De Mayo Universal Distractor®



- ◆ *New Removable "transformable"* carriage facilitates cleaning.



- ◆ Varus tilt control allows varying degrees of lateral movement by pushing the slide bar in increments from open to close.



Sterile Gel-infused
memory foam *imp*
Boot Pad &
Cohesive Wrap



Ordering Information

Description	Cat. No.
De Mayo V²™ Knee Positioner® System	114-KIT
INCLUDES: De Mayo V ² ™ / V ² E™ Knee Positioner® 30" Baseplate w/ Removable Carriage & Single Bar Lock Base	114-BP
De Mayo V ² ™ / V ² E™ Knee Positioner® Clamp	114-CL
Aluminum Boot w/ Distractor Block	803-ABD
De Mayo V²™ Knee Positioner® System	114-S KIT
INCLUDES: De Mayo V ² ™ / V ² E™ Knee Positioner® 25" Baseplate w/ Removable Carriage & Single Bar Lock Base	114-S BP
De Mayo V ² ™ / V ² E™ Knee Positioner® Clamp	114-CL
Aluminum Boot w/ Distractor Block	803-ABD
De Mayo V²™ Knee Positioner® System	114-907 KIT
INCLUDES: De Mayo V ² ™ / V ² E™ Knee Positioner® 30" Baseplate w/ Removable Carriage & Single Bar Lock Base	114-BP
De Mayo V ² ™ / V ² E™ Knee Positioner® Clamp	114-CL
Aluminum Boot w/ Distractor Block	803-ABD
De Mayo Universal Distractor®	907
De Mayo V²™ Knee Positioner® System	114-S-907 KIT
INCLUDES: De Mayo V ² ™ / V ² E™ Knee Positioner® 25" Baseplate w/ Removable Carriage & Single Bar Lock Base	114-S BP
De Mayo V ² ™ / V ² E™ Knee Positioner® Clamp	114-CL
Aluminum Boot w/ Distractor Block	803-ABD
De Mayo Universal Distractor®	907



De Mayo V²™ Knee Positioner is a trademark of Innovative Medical Products, Inc.
U.S. Patent No. 8,132,278

V²™ / V² E™ Clamp is a trademark of Innovative Medical Products, Inc.
U.S. Patent No. 7,003,827

De Mayo Universal Distractor® U.S. Patent No. 7,380,299

Proprietary Gel Foam Boot Pad property of Innovative Medical Products, Inc.

© July 2016 Innovative Medical Products, Inc.

ALL RIGHTS RESERVED

imp® **Innovative**
Medical Products, Inc.

87 Spring Lane, PO Box 8028

Plainville, CT 06062

PH: (860) 793-0391 FAX: (866) 459-1805

Toll Free PH: 1-800-467-4944

www.impmedical.com