

# Knee Positioning Triangles™

2 Triangles = 12 Positions

Simple, stable positioning  
for lower extremity surgery



Tibial Nail Position



Femoral Retrograde  
Nail Position



Lower Extremity  
Position



- ◆ Saves O.R. time.  
No fracture table setup.
- ◆ Radiolucent.
- ◆ Eases muscle tension during  
fracture reduction.
- ◆ Latex-free, sterile, single use  
Patient Protective Pads.
- ◆ Autoclavable.

## Ordering Information

Description	Cat. No.
Knee Triangle - Large (15x10x7)	193
Knee Triangle - Standard (12x7x7)	293
Knee Triangle Kit (1) 193, (1) 293	393 Kit
Sterile Patient Protective Pads 10/case	193-P

**imp® Innovative**  
Medical Products, Inc.  
A GLOBAL LEADER IN PATIENT POSITIONING

87 Spring Lane, PO Box 8028  
Plainville, CT 06062  
PH: (860) 793-0391 FAX: (866) 459-1805  
Toll Free PH: 1-800-467-4944  
[www.impmedical.com](http://www.impmedical.com)



Knee Positioning Triangles™ is a trademark of Innovative Medical Products, Inc.  
© 8/2013 Innovative Medical Products, Inc.  
ALL RIGHTS RESERVED